CLEARPOINT

NEURO

Flexible Adhesive Equipment Cover, 16x24

INSTRUCTIONS FOR USE

Table of Contents

1. Intended Use
2. Device Description
3. General Warnings and Precautions
4. Use Instructions
 - 4.1. Pre procedure
 - 4.2. Procedure - Use
 - 4.3. Post procedure

1. Intended Use

The ClearPoint Flexible Adhesive Equipment Cover is intended to protect equipment from contamination during procedures.

2. Device Description

Package Contents

NGS-EC-16X24-06-BR	Flexible Adhesive Equipment Cover, 16x24, 3 mil
20 Each	NON-STERILE

3. General Warnings and Precautions

Warning: Do not use any ferrous tools to trim or punch holes in the Equipment Cover while in the scanner room to prevent the danger of instruments flying into the scanner and possibly injuring people or equipment.

MRI Environment Compatibility:

This equipment cover material is MR Safe. It can be scanned safely under a static magnetic field.

4. Use Instructions

NOTE: The following instructions specifically address the method for covering Magnetic Resonance Imaging (MRI) Flexible Coils. Similar methods can be used to cover other instruments as necessary.

4.1. Pre procedure

CAUTION: Only non-magnetic materials may be used near the scanner.

- 4.1.1. When scissors or other cutting instruments are used in preparation of the Wrap, do not perform in the scanner room.
- 4.1.2. Clear an area approximately 2 ft x 2 ft.
- 4.1.3. Un-roll one sheet (16" x 24") and lay it sticky side up.

Note: There is a label on the 4th to last sheet as a reminder to order more Wrap.

Locating on Large Flex Coils

- 4.1.4. Center the flex coil over the Equipment Cover and press down gently on the coil and around the connector / cable area.
- 4.1.5. Place another piece of Equipment Cover sticky side down and press gently down onto the coil.

Locating on Small Flex Coils

- 4.1.6. Center the flex coil over one side of the Equipment Cover and press down gently on the coil and around the connector / cable area.
- 4.1.7. Fold the other side over on top of the coil, matching the perimeter, and press gently down onto the coil.

4.1.8. Perimeter attachment procedure.

- 4.1.8.1. Press the excess Equipment Cover around the perimeter firmly together.
- 4.1.8.2. Press carefully around the exiting cable to ensure a good seal.
- 4.1.8.3. Trim off an excess that is not desirable with scissors.

4.1.9. Coil gap locations procedure.

- 4.1.9.1. Gaps that will be pierced to allow fixation screws should have both sheets firmly pressed together.
- 4.1.9.2. Gaps that will not be breached should have only light pressure applied. The adhesion of the Wrap to the coil should be sufficient to hold it in place without adhering the Wrap to itself at every gap.

4.2. Procedure use

- 4.2.1. Use the flex coil as intended.
- 4.2.2. Create any openings that are necessary.
- 4.2.3. Re-seal any separation that might occur around the perimeter of the coil or in the gaps in the coil when noted.

4.3. Post Procedure

NOTE: Be careful when cutting around the coil to prevent cutting into the coil or any cabling coming out of the coil.

NOTE: Secure any cabling when removing Wrap in that area to prevent damage to the cable connection.

- 4.3.1. Remove the flex coil from the scanner room
- 4.3.2. Trim around the exterior to release the Equipment Cover where it has adhered to itself.
- 4.3.3. Gently pull the remainder of the Equipment Cover off the coil.
- 4.3.4. Dispose of the Wrap with biological waste per institution procedure.
- 4.3.5. Place the tube with remaining Equipment Covers back into the plastic bag / sleeve. Tuck the sleeve into the end of the tube to keep the wrap clean. The tube can be stored in this manner or in the carton it came in.

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ClearPoint Neuro, Inc